

Transcriptomics at the Genetics Core, Edinburgh Clinical Research Facility

Angie Fawkes & Katarzyna Hafezi

The **Genetics Core** provides support for transcriptomic projects through **RNA extraction**, **Library Preparation** and **Sequencing**. We can support many different protocols depending on your experimental needs. Below are some of our more popular workflows.

RNA Extraction

Blood – Paxgene (preferred) or Tempus
Tissue or Cells – Qiagen RNeasy Mini Kit
FFPE – Covaris truXTRAC FFPE Total NA Plus Kit
miRNA - Qiagen miRNeasy Mini Kit

Sample QC by Qubit & Agilent Fragment Analyzer, or Agilent Bioanalyzer, to determine RNA yield, RIN & DNA yield (to identify possible DNA contamination)

Library Preparation

Total RNA-seq (by ribodepletion)

Good – Information on coding & non-coding RNA; Can work with poor quality samples, including FFPE
Bad – Most expensive as requires a lot of reads
Can require 50-100M reads/sample

Input: 5ng-1µg & RIN: Any
 NEB NEBNext Ultra II Directional RNA Library Prep Kit for Illumina (rRNA depletion workflow)

or

Input: 250pg-10µg & RIN: Any
 Takara SMARTer Stranded Total RNA-Seq Kit v2 - Pico Input Mammalian

mRNA-seq (by polyA pull-down)

Good – Focus on coding regions
Bad – Requires good quality RNA
Can require 30-50M reads/sample

Input: 10ng-1µg & RIN: ≥7
 NEB NEBNext Ultra II Directional RNA Library Prep Kit for Illumina (polyA mRNA workflow)

or

Input: 2pg-200g & RIN: ≥7
 NEB NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina

Gene Expression Profiling

Good – Very affordable; works on poor quality samples; simple analysis
Bad – Information on level of transcripts only; no sequence data across the transcripts; less useful for isoforms.
Can require 5-10M reads/sample

Input: 0.5ng-500ng & RIN: Any
 Lexogen QuantSeq 3' mRNA-Seq Library Prep

miRNA-seq

A specific kit and sequencing run is required; miRNA data is not present in a standard RNA-seq or mRNA-seq run.
Can require 5-10M reads/sample

Input: ≥200ng & RIN: Any
 Perkin Elmer NEXTFLEX Small RNA-Seq Kit v3 for Illumina

or

Input: ≥1ng & RIN: Any
 Qiagen QIAseq miRNA Library Kit

Sequencing Platforms

Selecting the best platform for your experiment

2x NextSeq2000 & 1x iSeq for fast & flexible sequencing
 Choose the read length and number of reads you need per sample and we will do the rest.

Coming Soon – direct RNA sequencing on the ONT PromethION

Cycles (actual bp)	iSeq 100		NextSeq 2000			
	iSeq 4M Reads	P1 100M Reads	P2 400M Reads	P3 1200M Reads	P4 1800M Reads	
50 (88bp)				60Gb	90Gb	
100 (138bp)		10Gb	40Gb	120Gb	180Gb	
200 (238bp)			80Gb	240Gb	360Gb	
300 (327bp)	1.2Gb	30Gb	120Gb	360Gb	540Gb	
600 (638bp)		60Gb	180Gb			

Contact: A.Fawkes@ed.ac.uk
www.edin.ac/genetics-core